

Pfu DNA Polymerase

9K-001-0009 (250U) / 9K-001-0010 (500U)

store at -20°C

- . Pfu DNA Polymerase (2.5U/μl)
- . 10x Pfu reaction Buffer
- . 10 mM dNTP Mix (each 10mM)
- . Band Sharpener

Protocol

The recommended protocol for reaction mixture and temperature cycle in PCR is as below. You may modify the amount of template and Pfu DNA polymerase, extension time, annealing temperature, and cycle number according to the product size, primer T_m value, and template type.

| Reaction mixture (for 50μl reaction) | | Temperature cycle | |
|---|---------|-------------------|----------|
| 10X Pfu Taq buffer | 5.0 μl | 95°C 2 min | x1 |
| 10mM dNTP mix | 1.0 μl | 95°C 20 sec |] x10~40 |
| primer 1(10pmole/μl) | 3.0 μl | AT* 40 sec | |
| primer 2(10pmole/μl) | 3.0 μl | 72°C 1 min/kb | |
| template | X μl | 72°C 5 min | x1 |
| Band Sharpener | 0~20 μl | | |
| Pfu (2.5U/μl) | 0.50 μl | | |
| D.W. to | 50.0 μl | | |

AT*: Annealing Temperature

Adjust to the lower T_m between two primers

$$AT = T_m - (4 \sim 6^\circ C) \quad T_m = 2^\circ C \times (A+T) + 4^\circ C \times (G+C)$$

*It might be required to optimize PCR conditions depending on target size, T_m of primers, template nature, extension time, annealing temperature, enzyme quantity and cycle(s) numbers.

Technical Information

A. Template (temperature cycle)

- . Animal genomic DNA
 - 50-200 ng (25-35 cycles)
 - 10-50 ng (30-40 cycles)
- . Bacterial genomic DNA
 - 10-50 ng (20-25 cycles)
 - 1-5 ng (30-35 cycles)
- . Plasmid and lamda DNA
 - 1-5 ng (20-30 cycles)

B. Pfu DNA Polymerase

- For amplification >3 kb from animal genomic DNA, the amount of Pfu should be increased to 2 or 2.5U.
- LR DNA Polymerase is recommended for long-range PCR reaction.

C. Band Sharpener *

Band Sharpener is not necessary for regular PCR conditions. In case fragment include high GC region or hard to amplify complex secondary structures, please add Band Sharpener to a final concentration of 0.5x~2x (5-20μl for 50μl reaction) to reaction mixture, as optimization is required (see protocol below).

D. Primer design

- Primer can be designed using a primer design software or manually.
- Avoid repeated sequence at 3' end.
- In case 3'-end is G+C rich, the end have to be A or T.
- In case 3' end is A+T rich, the end have to be G or C.
- It is recommended that T_m of the designed primers is >64°C and AT >58°C.

E. Extension time

- Generally, extension should be performed at 2 min/kb.
- For amplification >5 kb, extension temperature should be assigned to 68°C.

Band Sharpener Optimization protocol

| Reaction mixture (Conc. of Band Sharpener) | Mix I (0.0 x) | Mix II (0.5x) | Mix III (1.0x) | Mix IV (1.5x) | Mix V (2.0 x) |
|---|------------------|------------------|-------------------|------------------|------------------|
| 10X Pfu buffer | 5.0 μl | 5.0 μl | 5.0 μl | 5.0 μl | 5.0 μl |
| 10mM dNTP mix | 1.0 μl | 1.0 μl | 1.0 μl | 1.0 μl | 1.0 μl |
| Primer 1 (10pmole/ μl) | 3.0 μl | 3.0 μl | 3.0 μl | 3.0 μl | 3.0 μl |
| Primer 2 (10pmole/ μl) | 3.0 μl | 3.0 μl | 3.0 μl | 3.0 μl | 3.0 μl |
| Template | X μl | X μl | X μl | X μl | X μl |
| Band Sharpener | 0 μl | 5.0 μl | 10.0 μl | 15.0 μl | 20.0 μl |
| Pfu (2.5U/μl) | 0.50 μl | 0.50 μl | 0.50 μl | 0.50 μl | 0.50 μl |
| Add D.W. to | 50.0 μl | 50.0 μl | 50.0 μl | 50.0 μl | 50.0 μl |

* **Note:** Band Sharpener included in the Kit is provided for PCR optimization purpose only. If your optimized PCR conditions include Band Sharpener, please order cat#: 9K-001-007.