



## **PRODUCT INFORMATION**

### **DNA EZ J1 (Endotoxin-Be-Gone) Solution**

**CAT.NO.:** DT71718

**Package:** 5ml

**Component:**

| Component                  | DT71718 |
|----------------------------|---------|
| Endotoxin-Be-Gone Solution | 5 ml    |
| Protocol                   | 1       |

**Storage**

Transportation at ambient temperature. Storage at 4°C. The period of validity see packaging.

**Features**

1. **Fast and simple. A single treatment takes 20 minutes.**
2. **High efficiency. 90% of the endotoxin can be removed by a single treatment; the endotoxin level is lower than 0.2 EU/ml after 3 treatments.**
3. **Easy to scale up. Flexible protocol.**

**Introduction**

Endotoxins, also called lipopolysaccharides (LPS), are major contaminants found in most regularly prepared biology samples such as plasmid, DNA, and proteins. The presence of small amounts of endotoxin in recombinant protein preparations can cause side effects in host organism such as endotoxin shock, tissue injury, and even death. Due to these reactions, it is essential to remove endotoxins from drugs, injectables, and other biological and pharmaceutical products.

The Solution Endotoxin-Be-Gone is an optimized solution for endotoxin removal from any biology samples in a solution form.

**Protocol**

1. **Add 50 µl of solution endotoxin-be-gone to 500 µl of DNA or protein solution.**
2. **Incubate on ice for 5 minutes.**
3. **Incubate at 37°C for 5 minutes.**
4. **Centrifuge at  $\geq 8000 \times g$  ( $\geq 10,000 \text{ rpm}$ ) for 5 minutes at room temperature.**
5. **Transfer the supernatant to a new endotoxin-free tube.**
6. **Repeat the procedure 2~3 times.**