

Product information

QF 24 V4
V1 Jan 2021

Rapid Animal Total RNA Extraction Kit

Catalog #: AT4181 / AT4182
Size: 50 preps / 250 preps
Storage: 4°C *

*: Check storage conditions for more details.

Product Description:

This kit uses a reagent similar to TRIzol for rapid RNA isolation from animal tissue, cells, and some of the plant tissue.

Features:

- Rapid and handy: the whole procedure takes approximately 10 minutes.
- Very low genomic DNA contamination (no detectable by RTPCR).
- Flexible, ideal for most animal tissue and cell culture.

Components:

Components	AT4181	AT4182
Solution A	50 ml	250 ml
Protocol	1	1

NOTE: Care must be taken when working with RNA. It is important to maintain an RNase-free environment starting with RNA sample preparation and continue through purification and analysis. Use RNase free tubes, tips, gels. Wear gloves at all times.

Procedures:

1. Preparation

- Estimate the amount of Solution A to be used. For 1 ml of Solution A, the maximum amount of tissue to be used is 100 mg; the maximum amount of cells to be used is 1×10^7 .
- Prepare enough chloroform, isopropanol, 75% ethanol and RNase-free water.

2. Homogenization

- For adherent cells*
Discard the culture medium, add 1 ml of Solution A to $\sim 10 \text{ cm}^2$ cultured cell, and mix gently by pipetting.
- For suspension cells*
Collect cells by centrifugation, discard the supernatant and add 1 ml of Solution A to $1-5 \times 10^6$ cells, and mix gently by pipetting. The amount of cells should not exceed 1×10^6 for fibroblasts or carcinoma cell.
- For fresh animal tissue*
Cut the tissue into pieces; add 1 ml of Solution A to 50-100 mg tissue, homogenate for 30 seconds using a homogenizer.
- For samples in RNAlocker*
Remove the liquid on the surface prior to cut, and then treat as fresh tissue.

3. Fraction Separation

- A. Add 0.2 volume of chloroform to the lysate, (eg. 0.2 ml chloroform for 1 ml lysate), vortex for 30 seconds.
- B. Centrifuge at no more than 12,000 x g (13,000-15,000 rpm) for 3 minutes at room temperature, deposit may form at the bottom of the centrifuge tube.
- C. Transfer the supernatant into a new 1.5 ml centrifuge tube.
- D. Repeat step 3 for 1-2 times for tissues rich in lipid.

4. Precipitation

- A. Add equal volume of isopropanol to the supernatant, vortex for 30 seconds.
- B. Centrifuge at no more than 12,000 x g (13,000-15,000 rpm) for 3 minutes at room temperature, deposit may form at the bottom of the centrifuge tube.
- C. Discard the supernatant carefully.

5. Primary Wash

- A. Add 1 ml of 75% ethanol to the RNA deposit, vortex for 30 seconds.
- B. Centrifuge at no more than 7,500 x g (13,000-15,000 rpm) for 3 minutes at room temperature.
- C. Discard the supernatant carefully.

6. Secondary Wash

- A. Repeat step 5 once.
- B. Short spin at no more than 7,500 x g (13,000 rpm), remove the residual ethanol carefully.

7. Dissolving the RNA

- A. Incubate the tube with lid open at room temperature for 1-2 minutes.
- B. Add 50-100 μ l of RNase-free water to dissolve the pellet, store RNA solution at -80°C

8. RNA integrity, yield and purity can be analyzed via electrophoresis and spectral analysis.

Storage:

Transport at room temperature, store all components at 4°C upon received.



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
NOT INTENDED FOR HUMAN OR ANIMAL USE.