

Technical information

SARS-CoV-2 RT-PCR Detection Kit


Catalog #: COV-2-RTPCR **ABI7500 Setup:** Page 1-2
BioRad CFX96 Setup: Page 2-4
Roche Lightcycler 480 II Setup: Page 4-8

Recommended Sample Type


The SARS-CoV-2 RT-PCR Detection Kit is used for the detection of SARS-CoV-2 viral RNA purified from various clinical samples. The recommended RNA samples are to be extracted from nasopharyngeal, or oropharyngeal swabs.

ABI7500 Setup


1. ABI 7500 Setup for Signal Collection:




2. No Passive Reference Dye Required:



3. Cycling Protocol:



4. Typical Results:



5. Additional Information:

Do we need to select Rox as a passive reference dye or should that be none?

– No passive reference required.

Do you have DMSO in your reaction mix?

– No.

Why does your kit recommend 45 cycles?

– Some positive samples show at Ct from 37-40, 45 cycles is for sufficient amplification.

6. Limitations of Use:

False negatives may be caused by:


- Unsuitable collection, handling and/or storage of samples.
- Sample outside of viraemic phase.
- Failure to follow handbook procedures.
- Use of unauthorised extraction kit or PCR platform.

False positives may be caused by:


- Unsuitable handling of samples containing high concentration of SARS-CoV-2 viral RNA or positive control template
- Unsuitable handling of amplified product.

BioRad CFX96 Setup


1. Open Bio-Rad CFX Manager Software and select “User-Defined” Run Type.




2. Under Run Setup > Protocol, select “Create New...” to create a new run protocol.




3. Create Cycling Protocol according to the datasheet. Click “OK”.




4. Under “Run Setup”, “Plate”, select “Create New...” to create a new Plate set up.




5. Use the “Select Fluorophores” Option to define channels: HEX (VIC), FAM and Cy5. Click “OK”.



6. Apply the three channels to all test wells by selecting the entire plate and checking “Load” for FAM, HEX and Cy5.




7. Under Run Setup > Start Run, Click on “Start Run” to start the run.




Roche Lightcycler 480 II Setup

1. Open Roche LightCycler 480 Software and select the correct “Plate Type” you are using.



2. Create a New Experiment.



3. Under “Run Protocol” > “Set Up”, select “4 Color Hydrolysis Probe” as the Detection Format. Then click on “Customize”. Check the channels for FAM, Red 610 (HEX) and Cy5. Click OK (✓).

The screenshot shows the 'Setup' window for a '4 Color Hydrolysis Probe' detection format. The 'Customize' button is highlighted. A dialog box titled 'Detection Formats' is open, showing the 'Active' filter combinations:

| Active | Filter Combination |
|-------------------------------------|-------------------------|
| <input type="checkbox"/> | Cyan 500 (440-488) |
| <input checked="" type="checkbox"/> | FAM (498-580) |
| <input checked="" type="checkbox"/> | Red 610 (533-610) |
| <input checked="" type="checkbox"/> | Cy 5 / Cy 5.5 (618-660) |

4. Setting up the run conditions. Use the (+) symbol to add 2 additional programs, and name each program accordingly. Enter the number of cycles required for each step. For “Analysis Mode”, make sure “Quantification” is chosen for the last Amplification step.

The screenshot shows the 'Programs' table with three programs:

| Program Name | Cycles | Analysis Mode |
|------------------------------------|--------|----------------|
| cDNA Synthesis | 1 | None |
| RT Inactivation & Pre-Denaturation | 1 | None |
| Amplification | 45 | Quantification |

The 'cDNA Synthesis Temperature Targets' table is also visible:

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 55 | None | 00:15:00 | 4.4 | 0 | 0 | 0 | 0 |

5. For each step, define the temperature and duration.

Click on “cDNA Synthesis” in the “Program Name” section, and then change the setting underneath according to the protocol.

LightCycler® 480 Software release 1.5.1.62 SP2

Instrument: LightCycler 480 / Standby (no MWP) Database: XDMS_R_new (Research) User: Applied Biological Materials

Window: New Experiment

Setup

Detection Format: 4 Color Hydrolysis Probe Block Size: 96 Plate ID: Reaction Volume: 20

Color Comp ID: Lot No: Test ID:

Programs

| Program Name | Cycles | Analysis Mode |
|------------------------------------|--------|----------------|
| cDNA Synthesis | 1 | None |
| RT Inactivation & Pre-Denaturation | 1 | None |
| Amplification | 45 | Quantification |

cDNA Synthesis Temperature Targets

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 55 | None | 00:15:00 | 4.4 | 0 | 0 | 0 | 0 |

Overview

Temperature (°C) vs. Estimated Time (h:mm:ss)

Apply Template End Program + 10 Cycles Start Run

6.1 Repeat step 5 for the other two programs.

LightCycler® 480 Software release 1.5.1.62 SP2

Instrument: LightCycler 480 / Standby (no MWP) Database: XDMS_R_new (Research) User: Applied Biological Materials

Window: New Experiment

Setup

Detection Format: 4 Color Hydrolysis Probe Block Size: 96 Plate ID: Reaction Volume: 20

Color Comp ID: Lot No: Test ID:

Programs

| Program Name | Cycles | Analysis Mode |
|------------------------------------|--------|----------------|
| cDNA Synthesis | 1 | None |
| RT Inactivation & Pre-Denaturation | 1 | None |
| Amplification | 45 | Quantification |

RT Inactivation & Pre-Denaturation Temperature Targets

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 95 | None | 00:00:30 | 4.4 | 0 | 0 | 0 | 0 |

Overview

Temperature (°C) vs. Estimated Time (h:mm:ss)

Apply Template End Program + 10 Cycles Start Run

6.2 For the “Amplification” Program, make sure that the “Acquisition Mode” is set as “Single”.

LightCycler@ 480 Software release 1.5.1.62 SP2

Instrument: LightCycler 480 / Standby (no MWP) Database: XDMS_R_new (Research)

Window: New Experiment User: Applied Biological Materials

Setup: Detection Format: 4 Color Hydrolysis Probe Block Size: 96 Plate ID: Reaction Volume: 20

Programs:

| Program Name | Cycles | Analysis Mode |
|------------------------------------|--------|----------------|
| cdNA Synthesis | 1 | None |
| RT Inactivation & Pre-Denaturation | 1 | None |
| Amplification | 45 | Quantification |

Amplification Temperature Targets:

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 95 | None | 00:00:15 | 4.4 | 0 | 0 | 0 | 0 |
| 60 | Single | 00:01:00 | 2.2 | 0 | 0 | 0 | 0 |

Overview: Temperature (°C) vs. Estimated Time (h:mm:ss) graph showing a 95°C hold followed by a 60°C amplification cycle.

7. Once the program is set up, you can enter sample names using the “Sample Editor” tab.

LightCycler@ 480 Software release 1.5.1.62 SP2

Instrument: LightCycler 480 / Standby (no MWP) Database: XDMS_R_new (Research)

Window: New Experiment User: Applied Biological Materials

Step 1: Select Workflow: Abs Quant, Rel Quant, Scanning, Color Comp, Tm, Melt Geno, Endpt Geno

Select Filter Combinations: 498-580, 533-610, 618-660

Abs Quant Units: []

Step 2: Select Samples

| Pos | Color | Repl Of | Sample Name |
|-----|-------------|---------|-------------|
| A1 | Blue | | Sample 1 |
| A2 | Red | | Sample 2 |
| A3 | Green | | Sample 3 |
| A4 | Magenta | | Sample 4 |
| A5 | Grey | | Sample 5 |
| A6 | Yellow | | Sample 6 |
| A7 | Brown | | Sample 7 |
| A8 | Cyan | | Sample 8 |
| A9 | Dark Green | | Sample 9 |
| A10 | Orange | | Sample 10 |
| A11 | Purple | | Sample 11 |
| A12 | Light Green | | Sample 12 |
| B1 | Blue | | Sample 13 |
| B2 | Red | | Sample 14 |
| B3 | Green | | Sample 15 |
| B4 | Magenta | | Sample 16 |
| B5 | Grey | | Sample 17 |
| B6 | Yellow | | Sample 18 |
| B7 | Brown | | Sample 19 |
| B8 | Cyan | | Sample 20 |
| B9 | Dark Green | | Sample 21 |
| B10 | Orange | | Sample 22 |
| B11 | Purple | | Sample 23 |
| B12 | Light Green | | Sample 24 |

Step 3: Edit Properties

Sample Name: []

Buttons: Apply Template, Configure Properties, Toggle View (Table), Reset All, Import, Export

8. To run the program, go back to the “Experiment” Tab and click “Start Run”.

The screenshot displays the LightCycler 480 Software interface. The 'Experiment' tab is selected and highlighted with a red box. The interface includes a top navigation bar with 'Run Protocol', 'Data', and 'Run Notes' tabs. Below this, there are fields for 'Detection Format' (4 Color Hydrolysis Probe), 'Block Size' (96), 'Plate ID', and 'Reaction Volume' (20). A 'Setup' section contains fields for 'Color Comp ID', 'Lot No', and 'Test ID'. The 'Programs' table lists the following steps:

| Program Name | Cycles | Analysis Mode |
|------------------------------------|--------|----------------|
| cDNA Synthesis | 1 | None |
| RT Inactivation & Pre-Denaturation | 1 | None |
| Amplification | 45 | Quantification |

The 'cDNA Synthesis Temperature Targets' table is also visible:

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 55 | None | 00:15:00 | 4.4 | 0 | 0 | 0 | 0 |

An 'Overview' graph shows the temperature profile over time. At the bottom, the 'Start Run' button is highlighted with a red box. Other buttons include 'Apply Template', 'End Program', and '+ 10 Cycles'.