

## Product information

QF 24 V4  
V1 Jan 2021

# Dual Source 75/150 Power Supply

**Catalog #:** EDK620  
**Size:** 1 UNIT  
**Storage:** 18-25°C\*

\*: Store unit at room temperature.

## Product Description:

The Dual Source 75/150 Power Supply is designed to power any combination of horizontal or vertical electrophoresis units. It generates constant voltage at settings of either 75 or 150 volts and has a maximum current of 300 milliamps. Safety features include a fused UL-approved, three-pronged power plug and a non-conductive casing. The unit is equipped with an easy-to-replace 1 Amp slow burn fuse.

## Specifications:

<b>Max Voltage:</b>	150 Volts
<b>Voltage Steps:</b>	75 or 150 Volts
<b>Max Current:</b>	300 Milliamps
<b>Type Output:</b>	Constant Voltage
<b>Lead Inputs:</b>	2 Set, Recessed, Color Coded
<b>Fuse:</b>	1.0 Amp
<b>Input Power:</b>	50 Hz, 110 Volts
<b>Connection:</b>	3-Wire Grounded Cord



## Operation:

Plug the unit into a grounded outlet and flip the switch to turn on the power. For the desired voltage, press the switch to the left or right (the left is 75 volts and the right is 150 volts). The lead inputs are color-coded. The black input is negative and the red input is positive.

## Maintenance:

This power supply includes a fuse tray that is located directly under the power cord socket at the rear of the unit (Fig 1). A flat head screw driver can be used to remove the tray (Fig 2) and replace the blown fuse when necessary (Fig 3).

The unit includes on spare fuse in the fuse tray. If this spare has already been used, the blown fuse can be replaced with a standard 5x20mm fuse.

