

BIO BASIC Worldwide



For more information on pricing, complete product line or to locate a Point of Sales near you, please visit our website or contact one of our Customer Service Representatives.

Email	✉	order@biobasic.com
Phone	☎	1 (905) 474-4493
Toll Free	☎	1 (800) 313-7224
Fax	☎	1 (905) 474-5794



pUCM-T Cloning Vector Kit

BS433 (1 µg)
BS434 (5 µg)
BS435 (20 Preps)
BS436 (100 Preps)
QF 24 TV4
CV1 2020

For Research Use Only

pUCM-T Cloning Vector Kit

Code: BS433 (1 µg)

BS434 (5 µg)

BS435 (20 Preps)

BS436 (100 Preps)

Description.....	1
Kit Contents.....	2
Storage.....	2
Preparation of Ligation Reaction.....	2
Ligation Reaction.....	3
Transformation.....	3-4
Screening.....	4
Identification of Transformants.....	4
Appendix.....	5-6

The T-Vector PCR Product Cloning Kit is suitable for cloning of PCR products with additional A at 3' end. The special ligation system provided by the kit enables customs to finish ligation in 1-2 hours.

pUCm-Tvector of our company is design for simplifying cloning of PCR products. Many thermal stable DNA polymerase, PCR products amplified by DNA polymerase such as Taq, Tth DNA produce additional A at 3' end which could be easily ligated to T vector with additional T.

pUCm-T of our company is a kind of novel pUC derivative T vector, the multiple restriction sites with most of them single site and adjusted beta-galactose reading frame make it easy for screening target clone through blue and white plaque. Specially designed two Pst I sites beside inserted fragments make it easy for screening target clone by Pst I digestion, at same time inserted fragments also could be screened by cheap and efficient restriction enzymes such as EcoR I and Hind III double digestion. Inserted fragments also could be sequenced using universal primers M13 and T7 promoter primer. In vitro transcription could be processed through site of T7 RNA polymerase promoter in pUCm-T.

Technical materials of our product are displayed at the end of protocol, the complete sequence of our T vector are same as pUC19 (GenBank Accession Number M77789) except difference at multiple cloning site.


B. Promoter and multiple cloning sites of pUCm-T vector

T7 Transcription Start

M13/pUC Sequencing Primer(-20) EcoR I Sac I Kpn I T7 Promoter Nde I

5'-G TTG TAA AAC GAC GGC CAG TGA ATT CGA GCT CGG TAC CGT AAT ACG ACT CAC TAT AGG GCG ACA TAT GAT

3'-C AAC ATT TTG CTG CCG GTC ACT TAA GCT CGA GCC ATG GCA TTA TGC TGA GTG ATA TCC CGC TGT ATA CTA



Xba I Sal I Pst I Pae I Hind III

TAG AGT CGA CCT GCA GGC ATG CAA GCT TGG CGT AAT CAT GGT CAT AGC TGT TTC CTG -3'

ATC TCA GCT GGA CGT CCG TAC GTT CGA ACC GCA TTA GTA CCA GTA TCG ACA AAG GAC -5'

Lac Z


M13/pUC Reverse Primer (-26)



PRODUCTS ARE FOR SCIENTIFIC RESEARCH ONLY
NOT INTENDED FOR HUMAN OR ANIMAL USE

Appendix

A. Map of pUCm-T vector



Kit Contents

Components	BS433	BS434	BS435	BS436
pUCm-T vector (50 ng/ul)	1 µg	5 µg	1 µg	5 µg
10x Ligation Buffer	-	-	50 µl	200 µl
50% PEG4000	-	-	50 µl	200 µl
T4 DNA ligase, 5U/ul	-	-	100 U	500 U
Sterilized ddH ₂ O	-	-	1 ml	1 ml

Storage

All contents of the kit should be kept in freezer after usage. The kit is stable for 9 months.

Preparation of Ligation Reaction

Purification of PCR products or not depends on quality of amplified product. If PCR products are very specific, purification of PCR products is not necessary. But if plasmids were used as template, it is necessary to purify PCR product, because plasmid template could form white colony after transformation. PCR products could be separated by agarose electrophoresis. PCR product gel extraction kit of our company (BS353/BS354) could recover >60 bp DNA fragment efficiently.

PCR products amplified by *Taq*, *Tth*, *AmpliTaq*, *KlenTaq* DNA polymerase bear additional A at 3' end. *Taq* DNA polymerase co-amplify with *Pfu*, *Pwo*, *Tli* or *Deep vent* DNA polymerase which possess 3'→5' exonuclease activity may bear an additional A at 3' end. PCR products possess additional A at 3' end could be ligated to pUCm-T. PCR products amplified by DNA polymerase with 3'→5' exonuclease activity is blunt end, cloning of this kind of fragments need to add additional A to blunt end.

A-Tailing Kit (BS513) of our company could be used for efficiently to add additional A to 3' end.

Ligation Reaction

1. In a standard ligation 10 µl reaction system, 50ng (1 µl) pUCm-T vector, 0.2 pmol PCR product, 1 µl 10 x ligation Buffer and 1 µl T4 DNA ligase, add ddH₂O to 10 µl final volume, usually it is not necessary to quantify PCR products. Ligation reaction could be manipulated as followings:

1µl	10x Ligation Buffer
1µl	50%PEG
1µl	pUCm-T vector
Xµl	purified PCR product
Yµl	Sterilized water
1µl	T4 DNA Ligase
<hr/>	
Final Volume	10µl

NOTE: Add T4 DNA ligase at last step.

2. Ligation for 1 hour to overnight at 16~23°C.

NOTE: Usually 1 hour is sufficient for normal experiment

Transformation

Suggest using Competent Cell Preparation Kit (BS523/524 or BS525/526)

1. Put 100 µl competent cell thawed on ice, gently mix well after complete thawed.
2. Add 5 µl ligation mixture, gently mix well. Put on ice for 30 minutes.
3. Heat shock at 42°C on water bath for 1 minute. Put on ice for 2 minutes.
4. Add 400 µl SOC culture shake at 37°C, 200-250 rpm for 1 hour.
5. Centrifuge at 4000 rpm for 5 minutes, discard 400 µl supernatant using pipette, resuspend the cells with rest of culture.

6. Plate bacteria at Amp⁺ plater which had been plated with 20 µl 100 mM IPTG and 100 µl 20 mg/ml X-gal.

NOTE: Quantity of bacteria used depends on efficiency of ligation and transformation

7. Discard surplus liquid after put the plate upward for 1 hour at 37 °C, put downward and culture overnight.

Screening

Screened transformants by blue/white colony.

When foreign DNA fragment were cloned into pUCm-T, reading frame of LacZ gene coding sequence was changed by insertion of foreign DNA, thus the activity of a-fragment was affected, recombinant clone on X-gal/IPTG plater display which colony at same time non recombinant display blue colony. Select which colony on IPTG/X-gal plater, transfer colony to Amp⁺ liquid with toothstick and culture overnight at 37°C.

Identification of Transformants

1. Extract plasmid from liquid culture of white colony, digest plasmid by PstI or other suitable restriction enzyme, and check the size of inserted fragment in order to decide whether the clone contains target fragment.
2. The sequence of target clone is determined by M13 universal primer or other suitable primer to further confirm the inserted fragment.