



## **PRODUCT INFORMATION**

### **DNase I –DD0099 (D0099)**

#### ***Product information for DD0099:***

##### **Description:**

Bio Basic Inc.'s DNase I (catalogue DD0099) is an endonuclease derived from bovine pancreas that will degrade double-stranded DNA in the presence of divalent cations, producing 3'-OH oligonucleotides.  $Mg^{2+}$  in the reaction solution causes the enzyme to produce nicks in double-stranded DNA, while in the presence of  $Mn^{2+}$ , DNase I cleaves both strands of the DNA.

DNase I is useful in nick translation for introducing single-stranded nicks that serve as primer sites for initiation of DNA synthesis and for cloning random DNA fragments by cleaving double-stranded DNA.

Bio Basic Inc.'s DNase I is a chromatographically pure preparation. It is offered as lyophilized powder with an approximately activity of 500 Kunitz U/mg. For greatest stability, it is important that DNase be dissolved at a concentration of at least 1mg/ml in 50% Glycerol with 20mM Tris-Cl, pH 7.5, and 1mM  $MgCl_2$ . This solution can be stored at -20°C for at least a year.

Recommended 1X Reaction Buffer: 40mM Tris-HCl (pH 8.0), 2.5mM (up to 10mM)  $MgSO_4$ , 1mM (up to 10mM)  $CaCl_2$

*Note: If starting material contains EGTA, EDTA, other chelating reagents, and/or high concentration of salts, higher concentration of Mg and Ca are recommended.*

Heat Inactivation: 10 minutes at 65°C in the presence of Stop Solution.

Inhibitors: EGTA; EDTA (7); salt concentrations >100mM will reduce DNase activity.

Molecular Weight: 31,000 Daltons.

Requirement:  $Ca^{2+}$  and  $Mg^{2+}$  or  $Mn^{2+}$

Source: Bovine pancreas.

Storage Temperature: Store at -20°C. Avoid exposure to frequent temperature changes. See the expiration date on the Product Information Label.

Stop Solution: 20mM EGTA (pH 8.0).

Unit Definition: One unit of DNase is defined as the amount required to completely degrade 1µg of lambda DNA in 10 minutes at 37°C in 50µl of a buffer containing 40mM Tris-HCl (pH 8.0), 2.5mM (up to 10mM)  $MgSO_4$ , 1mM (up to 10mM)  $CaCl_2$ . Under these assay conditions one unit of DNase activity is approximately equal to one Kunitz unit. See the unit concentration on the Product Information Label.

**Procedure:**

1. Add to an RNase-free PCR tube:

1 µg of RNA sample

49µl of 1X Reaction Buffer: 40mM Tris-HCl (pH 8.0), 2.5mM (up to 10mM) MgSO<sub>4</sub>, 1mM (up to 10mM) CaCl<sub>2</sub>

1 µl of DNase I, 1 unit/ml\*

\*Refer to the Certificate of Analysis for the lot specific activity. To dissolve DNase I at a concentration of at least 1unit/ml, storage buffer is recommended: 50% Glycerol, 20mM Tris-Cl, pH 7.5, and 1 mM MgCl<sub>2</sub>. This solution can be stored at -20°C for at least one year.

2. Incubate for 10~15 minutes at 37°C

3. To stop the reaction, add 1 µl of Stop Solution to bind calcium and magnesium ions and to inactivate the DNase I.

**Note:**


The Stop Solution (20 mM EGTA) must be added before heating to prevent metal (Mg/Ca) ion catalyzed hydrolysis of the RNA. Heat at 70 °C for 10 minutes to denature both the DNase I and the RNA.

This product should not be used for digestions longer than 15 minutes or for digestions at temperatures higher than 37°C, or the residual contaminating RNase activity will begin to degrade the RNA.

**Usage Notes:**

1. This DNase solution does not contain an RNase inhibitor.

2. Under different buffer conditions the amount of DNase required to completely digest a given amount of DNA may need to be empirically determined. For example, salt concentrations >100mM will reduce DNase activity.

**QC Results:**

Lane 1: Marker  
Lane 2: Negative control (no DNase I)  
Lane 3: DNA + 1 unit DNase I  
Lane 4: DNA + 2 units DNase I  
Lane 5: DNA + 3 units DNase I

lane 1 lane 2 lane3 lane4 lane 5