



Taq PCR Master Mix (2x, blue dye)

Kit Contents

Components	BS9295 1 ml	BS9296 5 x 1 ml
Taq PCR Master Mix (2x, blue dye)	1 ml	5 x 1 ml
Sterilized ddH ₂ O	1 ml	5 x 1 ml
MgCl ₂ (25 mM)	0.5 ml	1 ml
Protocol	1	1

Storage and stability

Transportation at 4°C. Store at -20°C.
Avoid repeated freeze-thaw cycles.

Introduction

Taq PCR Master Mix (2X, blue dye) is a ready-to-use solution containing Taq DNA polymerase, dNTP, MgCl₂, PCR buffer, PCR stabilizers, gel loading reagent and dye. PCR products can be loaded onto agarose gel directly. Optimized Taq PCR Master Mix (2X, blue dye) can amplify targets up to 5 kb in length from lambda DNA. Users only need to add a template, water and primers to set up a PCR reaction.

The Master Mix provides a final concentration of 1.5 mM MgCl₂ in the reaction mix. If a higher Mg²⁺ concentration is desired, please adjust Mg²⁺ concentration accordingly using 25 mM Mg²⁺ solution provided with the Master Mix.

Protocol

- Use the following guidelines to optimize the template concentration:
 - gDNA 1-10 µg/ml
 - plasmid DNA 0.1-1 µg/ml
- Optimizing the Primer Concentration: Primer concentrations in the range of 0.2–0.5 µM work for most PCR amplifications.
- Thaw reagents on ice.
- Prepare a reaction Master Mix using the following protocol:

Taq PCR Master Mix	25 µl	(1x)
DNA template	1 µl	(0.1-10 ng)
Primer F (10 µM)	2 µl	(0.4 µM)
Primer R (10 µM)	2 µl	(0.4 µM)
Nuclease-free ddH ₂ O	20 µl	
Total Volume	50 µl	

- Place the tube in a centrifuge and spin for 30–60 sec.
- Overlay PCR mixture with mineral oil when using a thermal cycler without a heated lid.
- Perform most PCR reactions using the following cycling program:

Pre-Duration	94°C	4 min	1 cycle
Duration	94°C	30 sec	30-35 Cycles
Anneal	45-68°C	30 sec	
Extend	72°C	1 kb/min	
Post-Extend	72°C	10 min	1 cycle
Hold	4 °C	N.A.	N.A.

Troubleshooting Guide

Possible Cause	Comment & Suggestions
<i>Little or no product.</i>	
Insufficient template DNA.	Increase template DNA to 0.1–1 µg/50 µl reaction.
Template DNA damaged or degraded.	Always use purified integrated high-quality DNA as templates.
Primer not match.	Review primer design.
Cycling program issue: Annealing temperature too low or too high; Extension time too short; Cycle number not enough.	Check the correct cycling conditions used. Adjust the cycling conditions to optimal. Perform a gradient PCR to determine the optimal annealing temperature if necessary.
Primer concentrations too low.	Use a primer final concentration of 0.2-0.5 µM.

Possible Cause	Comment & Suggestions
Mg ²⁺ concentration not optimal.	Optimize and adjust the concentration of Mg ²⁺ .
<i>Product is multi-banded.</i>	
Primer design not optimal.	Review primers and check them for specificity.
Primer degradation.	Check the concentration and quality of primer solutions.
Annealing temperature too low.	Increase the annealing temperature. Perform a gradient PCR to find the best annealing temperature.
<i>Product is smeared.</i>	
Too much starting template.	Reduce the amount of template DNA.
Mg ²⁺ concentration not optimal.	Optimize and adjust the concentration of Mg ²⁺ .
PCR of long	Increase Mg ²⁺

Possible Cause	Comment & Suggestions
fragments from large genomic DNA.	concentration. If the length of the product is 5 kb, please choose long PCR DNA polymerase.
Carry-over contamination.	Purify DNA templates. Preventing carry-over contamination with uracil-DNA glycosylase.



Email  order@biobasic.com
 Phone  1 (905) 474-4493
 T. Free  1 (800) 313-7224
 Fax  1 (905) 474-5794
 Web  www.biobasic.com