

Product information

QF 24 V4
V1 June 2020

dT-Tailing Kit

Catalog #: BS71217
Size: 100 preps
Storage: -20°C*

*: Product will be shipped with ice pack. Check storage conditions.

Product Description:

This kit is designed to add a single deoxythymidine (dT) residue at the 3' end of any blunt-end DNA fragments (normally a vector) based on the non-template-dependent terminal transferase activity of Taq polymerase. This T residue on the vector will hybridize to the single A overhang on the Polymerase Chain Reaction (PCR) product and increase the efficiency of ligation.

Storage:

Transportation on ice pack. Store at -20°C for 12 months.

Composition:

2 × dT-Tailing Buffer	2 x 1.25 ml
Taq DNA Polymerase (5 U/μl)	2 x 50 μl
Protocol	1

Procedures:

1. Blunt-end DNA fragments preparation

Blunt-end DNA fragments can be generated by digestion of plasmid (such as pUC19) with a blunt-end restriction enzyme (such as EcoR V) or by PCR using a DNA polymerase with proofreading activity (such as pfu). DNA should be purified by gel extraction, phenol:chloroform method or using a DNA clean up kit (Cat#BS367).

2. dT-Tailing reaction

In a thin-walled tube, add the following components:

Purified DNA fragments	0.2-2 μg
2 × dT-Tailing Buffer	25 μl
Taq DNA Polymerase (5 U/μl)	1 μl
ddH ₂ O	up to 50 μl

Mix thoroughly and incubate at 72°C for 2 hours.

3. Post Treatment

The dT-tailed DNA should be purified before further application, phenol:chloroform method is recommended.

- Add 50 μl of Phenol:Chloroform:Isoamyl Alcohol (25:24:1) to the reaction mixture and mix well.
- Centrifuge in a microcentrifuge at 12,000 rpm for 10 minutes to separate the phases.

- c) Transfer the aqueous phase (upper phase) to a fresh microcentrifuge tube.
- d) Add 50 μ l of Chloroform and mix well.
- e) Centrifuge at 12,000 rpm for 10 minutes to separate the phases.
- f) Transfer the aqueous phase (upper phase) to a fresh microcentrifuge tube.
- g) Add 100 μ l of ice-cold ethanol and mix well.
- h) Centrifuge at 12,000 rpm for 10 minutes to pellet the DNA.
- i) Discard the supernatant, add 100 μ l of ice-cold 70% ethanol and mix well. Centrifuge at 12,000 rpm to pellet the DNA.
- j) Repeat this step once.
- k) Discard the supernatant and air dry the DNA pellet.
- l) Dissolve the DNA pellet in ddH₂O or TE buffer.
- m) For ligation reactions, use approximately 25 to 50 ng of T-Vector.



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
NOT INTENDED FOR HUMAN OR ANIMAL USE.