

Product information

QF 24 V4
 V1 Jan 2020

Rapid Yeast Genomic DNA Extraction Kit

Catalog #: BS8227
Size: 50 preps
Storage: Mixed components storage*

*: Product will be shipped with ice pack. Check storage conditions.

Product Description:

The kit is designed for rapid small-scale extraction of high quality genomic DNA from yeast. Yeast cell wall is digested by Snailase (or lyticase, zymolyase). Whole cell is lysed by a special buffer and DNA is then precipitated and washed by alcohol. Purified DNA can be used for many downstream applications such as PCR, restriction enzyme digestion, hybridization and other applications.

Features:

- Rapid and Simple.
- High Quality of DNA. OD_{260}/OD_{280} of purified DNA is generally 1.8~1.9.
- No Toxic Substance. The kit does not contain toxic reagents.
- Easy to Scale Up.

Storage:

Upon receipt. Store kit at 4°C. For expiry date, please refer to kit label. Snailase should be stored at -20°C.

*Dilute 600 mg Snailase in 5 ml Snailase Storage Buffer before use. This is a Snailase Working Stock. Store Snailase Working Stock at -20°C.

Composition:

Universal Digestion Buffer	25 ml
Buffer PY	12 ml
TE Buffer	10 ml
Snailase Reaction Buffer	75 ml
Snailase Storage Buffer*	5 ml
Snailase	600 mg
Protocol	1

Procedures:

1. Collect 1.0 ml yeast culture ($\sim 1 \times 10^7$ cell) in a 1.5 ml Eppendorf tube and centrifuge at 10,000 x g (12,000 rpm) for 30 seconds. Discard supernatant completely.
2. Removal of yeast cell wall:
 - a) Enzymatic Digestion: Add 600 μ l Snailase Reaction Buffer, 1.2 μ l mercaptoethanol (not supplied in the kit) and 50 μ l Snailase Working Stock per 20 mg wet weight yeast in a 1.5 ml tube. Incubate at 37°C for 3 hours. Invert the tube periodically. If lyticase is used, add 50 μ l lyticase enzymatic storage buffer containing 300U or more lyticase per 20 mg wet weight yeast. Centrifuge at 3,000 x g (5,000 rpm) for 10 minutes. Discard the supernatant.
3. Add 400 μ l Universal Buffer Digestion, incubate at 65°C for 1 hour.

NOTE: To obtain RNA-free DNA, add 20 μ l Rnase A solution (20 mg/ml, not supplied in the kit) to the tube, mix thoroughly and incubate at room temperature for 5 minutes after 65°C incubation.
4. Add 200 μ l Buffer PY, mix by inverting, and incubate at -20°C for 5 minutes.
5. Centrifuge at 12,000 x g for 5 minutes at room temperature. Transfer the supernatant into a new 1.5 ml tube.
6. (Optional) Add 0.2 ml of chloroform to the supernatant, mix well by inverting 10 times. Centrifuge at 12,000 x g for 2 minutes. Carefully transfer the supernatant to a clean 1.5 ml tube.



BIO BASIC
Quality-Affordable Research

- 7.** Add equal volume of isopropanol (approx 0.3~0.5 ml) to the solution, mix well by inverting 5 times. Incubate at room temperature for 2~5 minutes. Centrifuge at 12,000 x g for 5 minutes, discard the supernatant carefully.
- 8.** Add 1 ml of pre-cooled 75% ethanol to the tube, mix well by inverting 10 times. Centrifuge at 12,000 x g for 1 minute, discard the supernatant.
- 9.** Repeat Step 8.
- 10.** Air-dry the pellet at room temperature with the lid open for 2~5 minutes.
- 11.** Add 50~200 μ l of TE buffer to dissolve DNA pellet. Keep at 4°C for a couple hours until DNA pellet is completely dissolved. Purified DNA is ready for use. Or keep at -20°C for long term storage.



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
NOT INTENDED FOR HUMAN OR ANIMAL USE.