

Product information

QF 24 V4
 V1 Feb 2020

TCA Precipitating Protein Kit

Catalog #: BSP011
Size: 50 preps
Storage: 2~-8°C*

*: Product will be shipped at ambient temperature. Check storage conditions. Products have one year expiration from time of purchase.

Product Description:

This kit is designed for precipitating and concentrating protein samples, even nanogram quantities, by TCA precipitation. The procedure is simple and rapid; it can be used to produce solubilized protein in diluents of your choice for SDS-PAGE electrophoresis or western blot. The kit contains (1) Trichloroacetic Acid 100% (6.1 N) and Deoxycholate Solution, (2) Wash Solution, (3) Dissolution buffer, and (4) Neutralization buffer. The kit is sufficient for 50 x 0.2 ml samples.

Storage and Transportation:

Transportation at room temperature. After receiving, store all reagents at 2-8°C. Always ensure that all caps are securely fastened to prevent volatilization.

Composition:

Precipitation Reagent A	2.5 ml
Wash Solution Buffer B	30 ml
Dissolution Buffer C	1 ml
Neutralization Buffer D	0.25 ml

Procedures:

1. Transfer 200 µl protein sample into a new 1.5 ml centrifuge tube.
2. Add 50 µl Precipitation Reagent A and vortex for 10 seconds.
3. Keep the mixture on ice or at 4°C for about one hour.
4. Centrifuge at 15000 rpm for 15 minutes at 4°C.
5. Carefully discard the supernatant and any residual liquid on inner surface of centrifuge tube. Keep the remaining pellet.
6. Add 600 µl cold Wash Solution (buffer B) and vortex for 10 seconds.
7. Keep the sample on ice or at 4°C for about 15 minutes, then centrifuge at 12000 rpm for 15 minutes at 4°C.
8. Carefully discard the supernatant and dry the pellets in fume hood or a SpeedVac.
9. Dissolve the pellets in 20 µl Dissolution Solution (buffer C) and vortex for 10 seconds. If the solution appears yellow, add 1-5 µl Neutralization solution (Buffer D) and vortex until the solution become blue.
10. Boil the sample for 5 minutes, and keep for SDS-PAGE electrophoresis.

Notes:

1. After step 5, it may be helpful to invert the centrifuge tube on filter paper to completely drain off residual liquid.
2. Always position microfuge tubes in the centrifuge with the same orientation, i.e. cap-hinge facing out-ward. This will allow the pellet to remain on the same side of the tube during repeated centrifugations and minimize the loss of protein. Additionally, always ensure the tip does not directly touch the mass of protein pellet when discarding the supernatant.
3. This kit can only be used for in vitro experiments.



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.