

Product information

QF 24 V4
V1 Feb 2025

Hotstart Taq DNA Polymerase

Catalog #: 3GHST81
Concentration: 5 U/ μ l
Size: 500U / 5x500U
Storage: -20°C

Product Description:

HotStart DNA Polymerase contains a proprietary antibody that blocks polymerase activity at low temperatures. During the initial denaturation step at 94°C, the antibody dissociates from DNA polymerase and restores enzyme activity. This feature significantly reduces non-specific product formations that would otherwise compete for reagent availability, thus offering higher specificity and improved yield of PCR products. PCR products amplified up to 6 kb in length with HotStart DNA Polymerase contain a single base (A) 3' overhang.

Source:

E.coli cells carrying a cloned pol gene from *Thermus aquaticus*.

Components:

Components	500U	5x500U(2,500U)
Hotstart Taq DNA Polymerase (5 U/ μ l)	100ul	500ul
5x Buffer	1.0ml	5x1.0ml

Storage Buffer Components*:

50 mM Tris-HCl (pH 8.0), 100 mM NaCl, 0.1 mM EDTA, 5 mM DTT, 50 % glycerol and 1.0 % Triton®X-100.

Unit Definition:

One unit of enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 mins at 70°C.

Shipping and Storage:

Upon arrival, HotStart DNA Polymerase should be stored at -20°C. Avoid repeated freeze-thaw cycles of all HotStart components to retain maximum performance.

*Intellectual property included

Basic PCR Protocol:

The following basic protocol serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (incubation times and temperatures, concentration of Taq Polymerase, primers, MgSO₄ and template DNA) vary and need to be optimized for each specific PCR.

All PCR experiments should be assembled in a nuclease-free environment. In addition, DNA sample preparation, reaction set-up and subsequent reaction(s) should be performed in separate areas to avoid cross-contamination.

A negative control reaction (omitting the template DNA) should always be performed in tandem with sample PCR to confirm the absence of DNA contamination.

1. Add the following components to a sterile 0.2ml PCR tube sitting on ice:

Components	Volume	Final Concentration
Template DNA	~100ng	~2ng
Forward Primer (10uM)	1 - 2.5ul	200 - 500nM
Reverse Primer (10uM)	1 - 2.5ul	200 - 500nM
5x Buffer	5ul	1X
dNTP Mixture (10 mM)	1ul	200 uM
Hotstart Taq Polymerase (5U/ul)	0.5ul	2.5 - 5U
Nuclease-free H ₂ O	25ul	-

- Mix contents of tube and centrifuge briefly.
- Incubate tube in a thermal cycler at 94°C for 10 mins to completely activate the HotStart DNA Polymerase and denature the template.
- Perform 30-35 cycles of PCR amplification as follows:
 - Denature:** 94°C for 30 sec
 - Anneal:** 45-72°C for 30 sec
 - Extend:** 72°C for 1 min/1kb template
- Incubate for an additional 5 mins at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.
- Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide or by using an Eco-DNA Dye (Cat #s DT81413, DT81414, DT81415, DT81417, DT81418). Use appropriate molecular weight standards.



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
NOT INTENDED FOR HUMAN OR ANIMAL USE.