

# BIO BASIC Worldwide



For more information on pricing, complete product line or to locate a Point of Sales near you, please visit our website or contact one of our Customer Service Representatives.

Email	✉	order@biobasic.com
Phone	☎	1 (905) 474-4493
Toll Free	☎	1 (800) 313-7224
Fax	☎	1 (905) 474-5794



One-Step M-MuLV RT-PCR Kit

BS665  
QF 24 TV4  
Cv2 2025

*For Research Use Only*

## One-Step M-MuLV RT-PCR Kit Code: BS665 (100 Preps)

Description.....	1
Features.....	2
Kit Contents.....	2
Storage.....	3
Quality Control.....	3
Procedure.....	3-4
Troubleshooting.....	5-6

### Description

Bio Basic's One-Step M-MuLV RT-PCR kit offers a unique system for performing RT-PCR in a single step, in a single tube. Traditionally, RT-PCR is performed in two reaction steps. In the initial reaction, primary-strand cDNA is reverse transcribed from total or poly(A) RNA. Then, in a separate reaction, the cDNA is amplified by PCR using a Taq DNA polymerase. This kit allows cDNA synthesis and PCR to be performed in a single optimized buffer, with both M-MuLV RT and Taq DNA polymerases. No additional reagents are required after the reaction is initiated. This method reduces the possibility of cross-contamination and provides a very convenient technique for gene expression.

### Low Specificity

#### A. Reaction conditions not optimal

- Optimize Mg concentration.
- Optimize primer.
- Optimize annealing temperature and extension time.
- Increase temperature of RT to 40-42°C.

#### B. Oligo-dT or random primers used for 1st strand synthesis

Use gene specific primers.

### Unexpected Bands

#### A. RNA contamination with genomic DNA

Pre-treat RNA with DNase I.



## Troubleshooting

### No Amplification Product

#### A. No cDNA Synthesis (temperature too low)

For cDNA synthesis step, incubate at 50°C.

#### B. RNA Contamination

Maintain RNase-free conditions, add RNase inhibitor.

#### C. Not enough starting material of RNA

Increase RNA concentration, use > 1ug of total RNA.

#### D. RNA has been damaged

Replace RNA.

#### E. RT inhibitors are present in RNA

Remove inhibitors in RNA preparation by additional 70% ethanol wash.

**NOTE:** Inhibitors include SDS, EDTA, spermidine, etc.

#### F. Annealing Temperature too high

Decrease temperature.

#### G. Extension time is too short

Set extension time for at least 1 minute per kb.

#### H. Cycle number is too low

Increase number of cycles.

## Features

- Rapid and convenient.
- Simple procedure. No transfer of sample between tubes is needed, which minimizes cross-contamination.
- Easy to detect more than one sample.
- 1.2 mM final concentration of magnesium in the reaction mix works well for most targets. If needed, the magnesium concentration can further be optimized (usually between 1.2-2 mM) with the 5 mM MgSO<sub>4</sub> solution provided.

**NOTE:** Care must be taken when working with RNA. It is important to maintain an RNase-free environment starting with RNA sample preparation and continuing through all purification and analysis. Use RNase free tubes, tips, and gels. Wear gloves at all times.

## Kit Contents

Components	BS665 (100 Preps)
10× One-Step RT-PCR Buffer	1 ml
M-MULV RT	50 µl
Taq DNA polymerase	50 µl
RNase Inhibitor (40 U/µl)	50 µl
RNase free ddH <sub>2</sub> O	2 ml
MgCl <sub>2</sub> (5 mmol/L)	1 ml
Protocol	1

## Storage

Transportation at 4°C. Store at -20°C. Valid for 1 year.

## Quality Control

This product is tested functionally using 12.5 pg of total HOX RNA as a template for the amplification of a 674-bp segment of beta-actin mRNA (40 cycles). A minimum of 25 ng of the RT-PCR product was obtained under these conditions.

## Procedure

1. Add RNase free ddH<sub>2</sub>O, RNA, 10X One Step RT-PCR Buffer, gene specific primers, M-MULV RT, Taq DNA polymerase into a sterile, nuclease-free tube on ice in the order indicated below:

Component	Volume (for 50 µl Reaction)	Final Concentration
10× One-Step RT-PCR Buffer	5 µl	1X
Forward primer (10 µmol/l)	2 µl	0.4 µmol/l
Reverse primer (10 µmol/l)	2 µl	0.4 µmol/l
M-MULV RT	0.5 µl	-----
Taq DNA Polymerase	0.5 µl	-----
RNase Inhibitor (40 U/µl)	0.5 µl	-----
Total RNA	x µl (<1 µg)	10pg – 1µg
RNase free water	Fill to 50 µl	-----

2. Mix gently, centrifuge briefly and place on ice.

3. Run the reactions using a thermocycler.

## Recommended Program Settings

42°C	10~30 minutes	} cDNA synthesis pre-denaturation
94°C	5 minutes	
94°C	30 seconds	} 30~40 cycles
40~65°C	30 seconds	
72°C	1 kb/minute	
72°C	10 minutes	

4. Setting the PCR parameters.
  - a. The optimal number of cycles depends on transcript abundance and template complexity, and must be determined empirically.
  - b. For experimental reactions, use 1 minute of extension time per kb.