



## Product information

# Taq Plus DNA polymerase

Cat. #	Size	Concentration
D0090	200U; 5x200U	5U/ul

### Product Description:

Taq Plus is a mixture of Taq and Proof Reading enzyme. Taq is a thermostable DNA polymerase isolated from a strain of *Thermus sp* (see product number D0089). Taq Plus is used to improve the reliability and yield of conventional primer extension reaction.

Taq Plus has **two following advantages:**

1. High fidelity with an error frequency  $1.6/10^6$  ( or  $0.0016/10^3$ ) during DNA synthesis.
2. Taq Plus increases the efficiency of polymerization reaction, resulting in a great percentage of extenuation reaction completion up to 10kb to 30kb.

Taq Plus has a temperature optimum between 72-78°C and remains > 95% active following 1-hour incubation at 95°C.

### 10xTaq Plus reaction buffer:

200mM TrisHCl (pH 8.8), 100mM KCl, 100mM (NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>, 20mM MgSO<sub>4</sub>, 1% Triton X-100, 1mg/ml bovine serum albumin (BSA).

**Note:** If not using this buffer, please note 0.1% TritonX-100 (final concentration) is a must to ensure high activity.

### Reaction Conditions:

**Note:** All reagents, including Taq Plus, should be mixed immediately before use.

DNA synthesis is performed in 100ul of mixture containing 20-200uM dNTPs, 0.3-1uM Primers, 0.1-0.250ng of template DNA, 10ul of 10 X reaction buffer and 2.5-5 units of Taq Plus. Mix the reaction gently, centrifuge briefly and then overlay with light mineral oil. Initially, denature the reaction by incubating at 95°C for 5 minutes and then cool to 40-68°C for 5 minutes to allow the primers to anneal to the template DNA.

### Optimization of DNA synthesis:

It is **important** to add the reaction components in the **following order:**

1. H<sub>2</sub>O
2. 10 x reaction buffer
3. dNTPs
4. DNA template and primers
5. Taq Plus



*Disclaimer: This product has not been licensed for use in the polymerase chain reaction (PCR) process for amplifying nucleic acids (US patent numbers 4683195 and 4683202 issued to Cetus).*



## Reaction Mixture Set Up:

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Add components, in the following order, into a thin-walled PCR tube. Keep all components on ice. The following control PCR reactions should be run in parallel to ensure that the Taq plus DNA polymerase is working properly.

Reagent	Final Concentration	Quantity	Reagent	Positive Control	Negative Control
Water (PCR--Grade)	---	variable	Water (PCR--Grade)	32.8ul	33.8ul
10x Taq Plus reaction buffer	1x	5ul	10x Taq Plus reaction buffer	5ul	5ul
2.5mM dNTP mixture	200uM of each	4ul	2.5mM dNTP mixture	4ul	4ul
Primer I, forward	0.1-1uM	variable	Primer I(10µM), forward	1ul	1ul
Primer II, reverse	0.1-1uM	variable	Primer II(10µM), reverse	1ul	1ul
Taq Plus DNA polymerase	1-1.5U/50ul	variable	Taq Plus DNA polymerase	0.2ul	0.2ul
Template DNA	See note 1	variable	Control DNA Template	1ul	---
Total Volume	---	50ul	Total Volume	50ul	50ul

3. Gently vortex the sample and briefly centrifuge to collect all drops from walls of the tube.
4. Overlay the sample with one-half of the total reaction volume of mineral oil or add an appropriate amount of wax.  
This step may be omitted if the thermo cycler is equipped with a heated lid.
5. Place samples in a thermo cycler and start PCR.

## Notes for the Components of the Reaction Mixture:

1. **Template DNA:** Usually the amount of template DNA is in the range of 0.01-1ng plasmid or phage DNA and 0.1-1ug for genomic DNA, for a total reaction mixture of 50ul.
2. **Primers:** The PCR primers are usually 15-30 nucleotides in length, longer primers provide higher specificity. The GC content of primer should be 40-60%. The primer should not be self-complementary or complementary to any other primer in the reaction mixture, and the melting temperature of flanking primers should not differ by more than 5°C. If the primer is shorter than 25 nucleotides, the approx. melting temperature (T<sub>m</sub>) is calculated using the formula as:  
 $T_m = 4(G+C) + 2(A+T)$ .
3. **dNTPs:** The final concentration of each dNTP in the reaction mixture is usually 200uM.
4. **Taq plus DNA polymerase:** Usually 1-1.5U of Taq plus DNA polymerase is used in the 50ul of reaction mix. Higher Taq plus DNA polymerase concentrations may cause synthesis of nonspecific products. However, if inhibitors are present in the reaction mix (e.g., if the template DNA used is not highly purified), higher amounts of Taq Plus DNA polymerase (2-3U) may be necessary to obtain a better yield of amplification products.
5. Usually the extending step is performed at 70-75°C. Taq Plus DNA polymerase exhibits lower extension rate than Taq DNA polymerase, so 2min extension time is recommended for every 1 kb to be amplified.
6. **Cycling conditions:** Usually denaturation for 0.5-2min at 94-95°C is sufficient; the optimal annealing temperature is 5C lower than the melting temperature of primer-template DNA duplex; Usually the extending step is performed at 70-75°C. Recommended extending time is 1min for the synthesis of PCR fragments up to 2kb. When larger DNA fragments are amplified, the extending time is usually increased by 1min for each 1kb.
7. **Number of cycles:** The number of PCR cycles depends on the amount of template DNA in the reaction mix and on the expected yield of the PCR product. For less than 10 copies of template DNA, 40 cycles should be performed. If the initial quantity of template DNA is higher, 25-35 cycles are usually sufficient.
8. **Final extending step:** After the last cycle, the samples are usually incubated at 72°C for 5-15min.