

## **Product information**

QF 24 V4 CV1 2018


## DNA Marker S Plus (100-5000bp) Ready-to-use

Catalog #:	SGM03 / SGM04
Size:	50 loadings / 5x50 loadings
Storage:	-20°C

## **Description:**

9 DNA fragments: 100, 250, 500, 750, 1000, 1500, 2000, 3000, 5000bp.

Usage:	6ul/load (containing 75ng of each DNA fragments and 150ng in 750bp band).	
<b>Recommended Loading:</b>	6ul (or 1ul/mm width of gel well).	
Storage:	-20°C	
Quality Control:	Agarose gel analysis shows that the bands between 100bp and 5000bp are accurate in size and distinguishable. The 750bp band is at least 2 times more intense than any other band.	



1.7% agarose gel

<mark>ک</mark> ۲

PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY. NOT INTENDED FOR HUMAN OR ANIMAL USE.