

## **Product Information**


Catalogue Number: GM303/GM304

Product Name: Low Range DNA Marker-A (25 bp-500 bp) Ready-to-Use

**Size:** 20 loadings (120ul)/ 5x20 loadings (5x120ul)

**Description: 9 DNA fragments: 25, 50, 75, 100, 150, 200, 300, 400, 500bp** The intensity of the 400bp band is increased to yield an internal reference indicator

Usage: 6ul /load containing 75ng of each DNA fragments and 150ng of 400bp band.



2.5% agarose gel

Storage: -20°C

Recommended Loading: 6µl (or 1µl/mm width of gel well)

Storage Buffer: 10mM Tris-HCl, pH 8.0, 1mM EDTA

## **Quality Control:**

- a) Agarose gel analysis shows that the bands are accurate in size and distinguishable.
- **b)** Incubation of  $6\mu l$  of the DNA KB-Ladder in restriction enzyme buffer overnight at  $37^{\circ}C$  shows absence of visible degradation.

QF 24 v2.0, 11/2011