

## 100bp DNA Ladder

9K-004-0002 (0.5mL)

store at -20°C



### Technical Information

The 100bp DNA ladder is suitable for sizing double - stranded DNA fragments. Ready to use without any handlings.

10 DNA fragments: 100, 200, 300, 400, 500, 600, 700, 800, 900 and 1000bp.

The intensity of the 500bp band is increased to yield an internal reference indicator.


Recommended loading: 6µL

Band Size Range: 100-1,000bp



#### Changes from 2014 Specifications:

1. Bands and intensities
2. Loading Quantity
3. % Agarose Gel



6µL loading on  
1.7% agarose gel