

PRODUCT INFORMATION

10-100µl Filter Tips

Description

10-100UL filter tip.

Pipette tips with filter inserts avoid contamination with liquids and aerosols during the pipetting. Filter tips minimize the danger of unwanted cross-contamination.

Features

- Filter tip
- 10-100UL
- Autoclavable



Item No.	Description	Package
FT224	Natural/filter/bag	1000/bag
FT224-RNS	Natural/filter/rack/sterile Certificated RNase, DNase-Free.	96/rack
FT224-RNSB	Natural/filter/rack/sterile Certificated RNase, DNase-Free.	96/rackX10

Technical Specifications

Item No.	FT224	FT224-NS	FT224-RNS
Package	1000/Bag	100/BagX10	96/RACK
Volume	10-100µl	10-100µl	10-100µl
Graduations	YES	YES	YES
Color	Natural	Natural	Natural
Plastic materials	Polypropylene	Polypropylene	Polypropylene
Filter	YES	YES	YES
Sterile	NO	YES	YES
DNase, RNase Free	NO	YES	YES
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C
Storage condition	Room temperature	Room temperature	Room temperature

Comparison of Related Tips



OF 24 Rev 1