




6x6 well Electrophoresis Apparatus

(Cat.#EDK626)



DNA electrophoresis for your whole class with just a single gel tank! Six groups of students can load their own individual gels and the six gels are run together in 30-40 minutes giving excellent results! Eliminate cumbersome gel tray taping and pour gels quickly and easily with our innovative gel tray sealing rubber end caps. The "6x6 well Electrophoresis Apparatus" is backed by a lifetime warranty and free technical support.

Features:

- New, sleek design.
- Contoured lid enhances visualization.
- Vented base.
- Color coded push tabs.
- Pour spout to pour out used buffer.
- Replaceable electrode modules.



For 6 Lab Groups



Cat. #EDK626

Product includes:

- One 6x6 well Electrophoresis Apparatus.
- Six 7x7cm Trays.
- Six 6-tooth Combs.
- Twelve Rubber End Caps.

