



Bio Basic Inc.

A world leader in serving science

## **Product Information**

### **Microtube Homogenizer**

#### ***Product information for BCM1201/BCM1202/BCM1203:***

#### **Introduction**

Homogenizer, ideal for use with microtubes. Capacity of up to 3 disposable screw-top microtubes. Ideal for samples that cannot be homogenized using a standard vortexer. 2 year warranty. For use with BCM1202 and BCM1203 prefilled tubes.



#### **Features:**

- Faster and more effective than tissue grinders
- Eliminates risk of cross-contamination that is present with hand-held homogenizers
- Powerful mixing action for lysis, grinding or homogenization
- Capacity of 3 tubes

#### **Ordering Information**

Product Code	Description
BCM1201	Microtube homogenizer, 115V
BCM1202	Triple-Pure starter kit, zirconium beads, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm bead tubes
BCM1203	Prefilled 5.0ml tubes, stainless steel beads, 2.8mm, acid-washed, 50pk

#### **Technical Data:**

<b>Speed</b>	2,800-4,000rpm	<b>Dimensions (WxDxH)</b>	6.9 x 8.3 x 5.3 inches 17 x 21 x 13.5 centimeters
<b>Timer</b>	3 sec to 3 mins	<b>Capacity</b>	3 x 2.0ml tubes (non-skirted)
<b>Weight</b>	5lbs/2.2kg	<b>Operating Temperature</b>	+4°C to +65°C
<b>Electrical</b>	115V, 50-60Hz, 50W	<b>Warranty</b>	2 Years

QF 24 Rev 0

V2 06/2012

20 Konrad Cres. Markham Ontario L3R 8T4 Canada  
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794  
Email: [order@biobasic.com](mailto:order@biobasic.com) Web: [www.biobasic.com](http://www.biobasic.com)