



Product Information Digital Water Baths

Product information for BCM1400/BCM1401/BCM1402/BCM1403:

Introduction

Digital water baths, with gabled lid. Gabled lid diverts condensation to front and back of bath. Lid is removable. BCM1400 has a flat, non-removable lid. 8L and 12L models have a drain for easy removal of water. All models have seamless, stainless steel interiors which are rust-resistant and have no corners to collect debris or contaminants. 2 year warranty.



Features:

- Hinged, removable lid (except for BCM1400 – lid cannot be removed)
- Seamless stainless steel chamber
- Digital controls

Ordering Information

Product Code	Description
BCM1400	2L digital water bath, 115V
BCM1401	4L digital water bath, 115V
BCM1402	8L digital water bath, 115V
BCM1403	12L digital water bath, 115V

Technical Data:

	BCM1400	BCM1401	BCM1402	BCM1403
Temperature Range	Ambient +5°C to 100°C			
Temperature Accuracy	±0.2°C at 37 °C		±0.5°C at 37 °C	
Temperature uniformity	±0.2°C at 37 °C			
Exterior Dimensions (WxDxH)	6.5 x 8 x 6 inches	7.7 x 13.9 x 9 inches	17.5 x 10.5 x 12.4 inches	18 x 10.5 x 16.2 inches
Chamber Dimensions (WxDxH)	5.5 x 6 x 4 inches	5.3 x 9.4 x 5.5 inches	12.8 x 10.4 x 5.9 inches	12.8 x 12.4 x 7.9 inches
Weight	4.5 lbs / 2kg	15 lbs / 6.5 kg	17 lbs / 8 kg	17 lbs / 8 kg
Electrical	115V, 50-60 Hz			
Warranty	2 years			

QF 24 Rev 0

V2 06/2012

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com