




# Product Information

## Tube Rollers



### Product information for BCM1608/BCM1609:

#### Introduction

Fixed speed tube rollers, for mixing contents of tubes. Compatible with tubes and bottles with diameters from 8mm to 200mm. The right side of each roller is raised and lowered during each revolution to provide a gentle up-down mixing motion in addition to the rolling mixing motion.

#### Features:

- Gentle rolling motion
- 5° tilting motion

#### Ordering Information

Product Code	Description
BCM1608	Tube roller with 5 rollers, 115V
BCM1609	Tube roller with 10 rollers, 115V

#### Technical Data:

<b>Speed</b>	38rpm±10%	<b>Dimensions (WxDxH)</b>	<b>BCM1608:</b> 18 x 7 x 3.5 inches 45 x 18 x 9 centimeters <b>BCM1609:</b> 18 x 14 x 3.5 inches 45 x 36 x 9 centimeters
<b>Motion</b>	Simultaneous roll & tilt (5°)	<b>Roller Spacing</b>	5mm, on 35mm centers
<b>Vessels Accepted</b>	8mm-200mm diameter	<b>Operating Temperature</b>	+4°C to +40°C
<b>Weight</b>	<b>BCM1608:</b> 11lbs/5kg <b>BCM1609:</b> 17.6lbs/8kg	<b>Warranty</b>	2 Years
<b>Electrical</b>	100V, 50-60Hz, 250W		

#### Capacity:

Vessel	BCM1608	BCM1609
12 x 75mm tubes	12	26
16 x 100mm tubes	8	18
16 x 125mm tubes	8	18
15ml screw cap tubes	8	18
50ml screw cap tubes	8	18
80cm <sup>2</sup> roller bottles	1	2