



Product Information Junior Orbital Shaker

Product information for BCM1622:

Introduction

Orbital shaker with flat mat. Ideal for aeration and mixing. Compact model. Speed can be adjusted. 2 year warranty.

Features:

- Adjustable speed
- Includes non-slip rubber mat platform
- Safe for incubator or cold room use
- Compatible with BCM1623 magnetic clamp platform for flasks and tube racks (pictured to right; not included with BCM1622)



Ordering Information

Product Code	Description
BCM1622	Shaker JR with rubber mat platform, 9.5" x 8", 100-240V, US plug

Technical Data:

Speed Range	50-300rpm; adjustable in increments of 10rpm	Dimensions (WxDxH)	11 x 10.7 x 5inches 28 x 27 x 13 centimeters
Orbit	19mm	Timer	1min to 99hr,59min
Load Capacity	4 x 1L, or 9lbs (4kgs)	Platform Size	9.5 x 8 inches 24 x 20 cm
Weight	12lbs/5.5kg	Operating Temperature	+4°C to +65°C
Electrical	Universal Voltage 100-240V, 50-60Hz, 250W	Warranty	2 Years