



## Product Information Orbital Shaker

### *Product information for BCM1624:*

#### Introduction

Orbital shaker with flat mat. Ideal for aeration and mixing. High Capacity. Speed can be adjusted. 2 year warranty.

#### Features:

- Adjustable speed
- Includes non-slip rubber mat platform
- Safe for incubator or cold room use
- Compatible with BCM1623 magnetic clamp platform for flasks and tube racks (pictured to right; not included with BCM1624)



#### Ordering Information

Product Code	Description
BCM1624	Shaker with rubber mat platform, 13" x 12", 115V

#### Technical Data:

<b>Speed Range</b>	30-300rpm; adjustable in increments of 10rpm	<b>Dimensions (WxDxH)</b>	11 x 13 x 3 inches 28 x 33 x 7.5 centimeters
<b>Orbit</b>	19mm	<b>Timer</b>	1min to 48hr/continuous
<b>Load Capacity</b>	5 x 1L, or 14lbs (6.4kgs)	<b>Platform Size</b>	13 x 12 inches 33 x 30 cm
<b>Weight</b>	25lbs/11.5kg	<b>Operating Temperature</b>	+4°C to +65°C
<b>Electrical</b>	Universal Voltage 100-240V, 50-60Hz, 250W	<b>Warranty</b>	2 Years