



Product Information Heated Microplate Vortexer

Product information for BCM1630/BCM1631:

Introduction

Temperature controlled microplate vortexer. Accepts a variety of microplate formats. 2 year warranty.

Features:

- Adjustable temperature, time and speed
- Compact footprint
- Heated lid prevents condensation



Ordering Information

Product Code	Description
BCM1630	Mixer, heated plate vortexer, 2 plate capacity, 115V
BCM1631	Mixer, heated plate vortexer, 4 plate capacity, 115V

Technical Data:

Speed Range	BCM1630: 100-1200rpm BCM1631: 100-1500rpm 10rpm increments	Dimensions (WxDxH)	BCM1630: 11 x 10.6 x 5.5 in. 28 x 27 x 14 cm	BCM1631: 12.6 x 13.8 x 7.3 in. 32 x 35 x 18.5 cm
Temperature Range	Ambient 5°C to +70°C, 0.1°C increment	Temperature Accuracy	±0.2°C	
Orbit	2mm	Temperature Uniformity	±0.5°C	
Weight	18lbs/8kg	Timer	1 min to 99hr, 99min / continuous	
Electrical	115V, 50-60Hz, 500W	Warranty	2 Years	