



# Product Information

## Vortexer with Heating and Cooling

Product information for  
BCM1632/BCM1633/BCM1633-MP:

### Introduction

Temperature controlled vortexer. Display shows both actual and selected values. 2 year warranty.

### Features:

- Adjustable temperature and speed
- Compact footprint
- Optional microplate block available
- BCM1632 has a sequence link function for linking or repeating programs



### Ordering Information

Product Code	Description
BCM1632	Shaker with heating and cooling, 115V
BCM1633	Shaker with heating, 115V
BCM1633-MP	Block for one microplate, for BCM1632 and BCM1633

### Technical Data:

<b>Speed Range</b>	200-1500rpm, 3mm orbit	<b>Dimensions (WxDxH)</b>	8.5 x 11 x 8.5 inches 21 x 29 x 21 centimeters
<b>Temperature Range</b>	<b>BCM1632:</b> -20°C to +100°C <b>BCM1633:</b> +5°C to +100°C	<b>Temperature Accuracy</b>	±1.0°C
<b>Heating Ramp Rate</b>	5.0°C/min	<b>BCM1632 Cooling Ramp Rate</b>	2.5°C/min
<b>Weight</b>	18lbs/8kg	<b>Warranty</b>	2 Years
<b>Electrical</b>	115V, 50-60Hz, 500W		