



Product Information

Large Platform Rocker

**Product information for
BCM1643/BCM1644/BCM1645:**


Introduction

Variable platform rocker with large work surface. Speed and tilt can be adjusted. Optional accessories are available. 5 year warranty.

Features:

- Adjustable tilt and speed
- Large work surface
- Optional stacking platform to increase working area
- Safe for incubator or cold room use

Pictured with optional stacking platform, BCM1644



Ordering Information

Product Code	Description
BCM1643	Rocker with flat mat, 115V
BCM1644	Stacking platform for BCM1643 with flat mat, 3.5" clearance
BCM1645	4 cord bungee pack, 7" to 18" working length

Technical Data:

Speed Range	3-80rpm	Dimensions (WxDxH)	16 x 13 x 6 inches 40.5 x 33 x 15.25 centimeters
Tilt Angle	0° to ±20°	Platform Size	16 x 12 inches 40.5 x 30 cm
Load Capacity	16lbs 7.5kg	Operating Temperature	0°C to +65°C
Weight	20lbs/9kg	Warranty	5 Years
Electrical	115V, 50-60Hz, 50W		