



Product Information

Mini Microcentrifuge

Product information for BCM1701/BCM1702/BCM1703/BCM1704/BCM1705/BCM1706:

Introduction

Compact, personal microcentrifuges ideal for quick spindowns of microtubes and PCR tubes. Includes optional PCR strip rotor, 0.2ml tube adapters and 0.5ml tube adapters. 1 year warranty.




Features:

- Click-n-lock rotor attachment system
- Eight position microtube rotor,
- Storage compartment for PCR rotor and 0.5ml Tube adapters
- Decelerates in one second

Bonus accessories, included:

- 0.2ml Tube Adapters
- 0.5ml Tube Adapters
- Four position 0.2ml PCR Strip Rotor
- Holds 32 x 0.2ml PCR tubes



Storage Compartment



Includes Both Rotors

Ordering Information

Product Code	Description
BCM1701	Mini Centrifuge, Blue Lid, with 2 rotors, 100-240V, US Plug
BCM1702	Mini Centrifuge, Clear Lid, with 2 rotors, 100-240V, US Plug
BCM1703	Mini Centrifuge, Green Lid, with 2 rotors, 100-240V, US Plug
BCM1704	Mini Centrifuge, Purple Lid, with 2 rotors, 100-240V, US Plug
BCM1705	Mini Centrifuge, Red Lid, with 2 rotors, 100-240V, US Plug
BCM1706	Optional rotor, 2 x 8 position PCR strip capacity

Technical Data:

Speed	6,000rpm/2,000xg	Dimensions (WxDxH)	4.5 x 5.9 x 4.5 inches 11.4 x 15 x 11.4 centimeters
Capacity	8 x 1.5-2ml Tubes 32 x 0.2ml PCR Tubes 4 x PCR Strips (8 x 0.2ml)	Weight	2.2lbs/1kg
Electrical	Universal Voltage 100-240V, 50-60Hz, 250W	Warranty	1 Year

QF 24 Rev 0



Bio Basic Inc.

A world leader in serving science

QF 24 Rev 0

V2 06/2012

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com