



Product Information **Mini Microcentrifuge**

Product information for BCM1707/BCM1708:

Introduction

Compact, personal microcentrifuges ideal for quick spindowns of microtubes and PCR tubes. Includes optional PCR strip rotor and 0.5ml tube adapters. 1 year warranty.

Features:

- Rotors can be changed in seconds
- Six position microtube rotor
- Decelerates in seconds

Bonus accessories, included:

- 0.5ml Tube Adapters
- Two position 0.2ml PCR Strip Rotor
- Holds 16 x 0.2ml PCR tubes



Ordering Information

| Product Code | Description |
|---------------------|---|
| BCM1707 | Red Mini Centrifuge with 2 rotors, 100-240V, US Plug |
| BCM1708 | Blue Mini Centrifuge with 2 rotors, 100-240V, US Plug |

Technical Data:

| | | | |
|-------------------|---|---------------------------|--|
| Speed | 6,600rpm/2,200xg | Dimensions (WxDxH) | 5 x 5.4 x 3.9 inches 12.6 x 13.8 x 10 centimeters |
| Capacity | 6 x 1.5-2ml Tubes 16 x 0.2ml PCR Tubes 2 x PCR Strips (8 x 0.2ml) | Warranty | 1 Year |
| Electrical | 100-240V, 50-60Hz, 0.1A | | |