



Product Information Mini Centrifuge

Product information for BCM1746/BCM1747/BCM1748:

Introduction

Low profile mini centrifuge, with nearly twice the capacity of standard mini centrifuges. Includes dual purpose rotor, compatible with both strips and tubes. 1 year warranty.

Features:

- Quiet motor
- Multi-use rotor, compatible with both strips and tubes
- Starts and stops with opening and closing of the lid



Ordering Information

Product Code	Description
BCM1746	12 Mini Centrifuge, 100-240V (US Plug)
BCM1747	Adapters for 0.5ml tubes, package of 6
BCM1748	Adapters for 0.2ml tubes, package of 6

Technical Data:

Speed	5,500rpm/2,000xg	Dimensions (WxDxH)	5.5 x 7.9 x 4.4 inches 14 x 20 x 11.2 centimeters
Capacity	12 x 1.5-2ml Tubes 32 x 0.2ml PCR Tubes 4 x PCR Strips (8 x 0.2ml)	Weight	2.2lbs/1kg
Electrical	Universal Voltage 100-240V, 50-60Hz, 300W	Warranty	1 Year