



## Technical Data Sheet

### Recombinant Murine KC/CXCL1 (rMuKC/CXCL1)

#### *Murine KC/CXCL1*

The three GRO cDNAs encode 107 amino acid precursor proteins from which the N-terminal 34 amino acid residues are cleaved to generate the mature GROs. There are no potential N-linked glycosylation sites in the amino acid sequences. GRO expression is inducible by serum or PDGF and/or by a variety of inflammatory mediators, such as IL-1 and TNF, in monocytes, fibroblasts, melanocytes and epithelial cells. In certain tumor cell lines, GRO is expressed constitutively. Similar to other alpha chemokines, the three GRO proteins are potent neutrophil attractants and activators. In addition, these chemokines are also active toward basophils. All three GROs can bind with high affinity to the IL-8 receptor type B.

Catalog Number:	RC332-12
Source:	<i>Escherichia coli</i> .
Molecular Weight:	7.8 kDa, a single, non-glycosylated polypeptide chain containing 72 amino acids.
Quantity:	1.0 mg
AA Sequence:	APIANELRCQ CLQTMAGIHL KNIQSLKVLP SGPHCTQTEV IATLKNAREA CLDPEAPLVQ KIVQKMLKGV PK
Purity:	>97% by SDS-PAGE and HPLC analyses.
Biological Activity:	Fully biologically active when compared to standard, determined by its ability to chemoattract total human neutrophils using a concentration range of 10.0-100.0 ng/ml.
Formulation:	Lyophilized from a 0.2µm filtered concentrated solution in 20mM PB, pH 7.4, 150mM NaCl.
Endotoxin:	Less than 1EU/µg of rMuKC/CXCL1 as determined by LAL method.
Reconstitution:	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml. Stock solutions should be apportioned into working aliquots and stored at ≤-20°C. Further dilutions should be made in appropriate buffered solutions.



Bio Basic Inc.

A world leader in serving science

---

Storage:

This lyophilized preparation is stable for several weeks at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Usage:

This material is offered by Bio Basic Inc for research, laboratory or further evaluation purposes. NOT FOR HUMAN USE.