



Product Information

EZ-LINE Cell Culture Dishes, Treated

Product information for SP31146, SP31171, SP31212:

Description

Cell culture dishes for tissue culture work. All dishes are produced with optically clear, high quality crystalline polymers. Growth surface is treated for optimal cell attachment and growth and designed for reduced cross-contamination.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- TC treated
- High strength



Item No.	Description	Package
SP31146	EZ-LINE 35x10mm Cell Culture Dish, TC Treated, sterile, 10 per pack, 500 per case, External Grip	500/case
SP31171	EZ-LINE 60x15mm Cell Culture Dish, TC Treated, sterile, 10 per pack, 500 per case, External Grip	500/case
SP31212	EZ-LINE 100x20mm Cell Culture Dish, TC Treated, Sterile, 10 per pack, 200 per case, External Grip	200/case

Technical Specifications

Item No.	SP31146	SP31171	SP31212
Quantity per case	500/case	500/case	200/case
Packaging	10 Pack	10 Pack	10 Pack
Volume	3mL	5mL	12.5mL
Growth Area	9.40cm ²	21.5cm ²	57.0cm ²
Color	Clear	Clear	Clear
Plastic materials	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes
TC Treated	Yes	Yes	Yes

Storage

Room temperature.

QF 24 Rev 0