



Product Information EZ-LINE Multi-Channel Reservoir Trays

Product information for SP33161, SP32113, SP32119, SP32123:

Description

Single and multi-channel reservoirs for cell culture work. Wells have notched volume labels and sloping sides for easy pipetting. Volumes refer to each individual channel.



Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- Notched volume labels
- High strength

Item No.	Description	Package
SP33161	EZ-LINE Reservoir, 50ml, 1 Channel, sterile, indiv. wrapped, 50 per case	50/case
SP32113	EZ-LINE Reservoir, 25ml, 2 Channel, sterile, indiv. wrapped, 50 per case	50/case
SP32119	EZ-LINE Reservoir, 7ml, 8 Channel, sterile, indiv. wrapped, 100 per case	100/case
SP32123	EZ-LINE Reservoir, 4ml, 12 Channel, Sterile, indiv. wrapped, 100 per case	100/case

Technical Specifications

Item No.	SP33161	SP32113	SP32119	SP32123
Quantity per case	50/case	50/case	100/case	100/case
Packaging	Individually wrapped	Individually wrapped	Individually wrapped	Individually wrapped
Volume/Channel	50mL	25mL	7mL	4mL
Color	White	White	White	White
Plastic materials	Polystyrene	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes	Yes

Storage

Room temperature.

QF 24 Rev 0