



Product Information

EZ-LINE T-Flask, with Vented Filter Cap, Treated

Product information for SP81136, SP81186, SP82286:

Description

T-flasks with vented filter caps for tissue culture work. All flasks are produced with optically clear, high quality crystalline polymers. Growth surface is treated for optimal cell attachment and growth and designed for reduced cross-contamination.



Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- TC treated
- High strength

Item No.	Description	Package
SP81136	EZ-LINE 25cm T-Flask, TC Treated, Vented Filter Cap sterile, 5 per pack, 200 per case	200/case
SP81186	EZ-LINE 75cm T-Flask, TC Treated, Vented Filter Cap sterile, 5 per pack, 100 per case	100/case
SP82286	EZ-LINE 24-well TC Plate, TC Treated, Vented Filter Cap Sterile, 5 per pack, 40 per case	40/case

Technical Specifications

Item No.	SP81136	SP81186	SP82286
Quantity per case	200/case	100/case	40/case
Packaging	5 Pack	5 Pack	5 Pack
Volume	7mL	25mL	50mL
Growth Area	25.0cm ²	75.0cm ²	175.0cm ²
Color	Clear, Green Cap	Clear, Green Cap	Clear, Green Cap
Plastic materials	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes
TC Treated	Yes	Yes	Yes

Storage

Room temperature.

QF 24 Rev 0