



Product Information

Micro UV Quartz Cell With Frosted Walls and With Lid

Product information for SQ244/SQ264/SQ274:

Description



Item No.	H x Lx W (mm)	Path length	Inside width	Volume	UV Range
SQ244	45x12.5x12.5	10mm	2mm	0.7ml	190nm-2500nm
SQ264	45x12.5x12.5	10mm	4mm	1.4ml	190nm-2500nm
SQ274	45x12.5x12.5	10mm	5mm	1.7ml	190nm-2500nm

Storage

Room temperature.