



# Product Information

## 50ml Skirted-Self-Standing Shaped With Screw Cap & Graduated

### Product information for TK943-TK944:

#### Description:

50ml Polypropylene Skirted-Self-Standing Shaped With Screw Cap & Graduated. Unlike polystyrene tube, Polypropylene tubes are chemical inert. Autoclavable. Tubes with screw cap are suitable for cell culture.

#### Features:

- Polypropylene tube
- 50ml
- screw cap
- Suitable for cell culture
- Autoclavable



Item No.	Description	Package
TK943-Y	Tube & yellow cap	25/bag
TK943-G	Tube & green cap,	25/bag
TK944-YS	Tube & yellow cap, sterile	25/bag
TK944-GS	Tube & green cap, sterile	25/bag

#### Storage

Room temperature.

#### Technical Specifications

	TK943-Y	TK943-G	TK943-YS	TK943-GS
Tube Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Tube Color	Natural	Natural	Natural	Natural
Capacity(ml)	50ml	50ml	50ml	50ml
Bottom Style	Conical	Conical	Conical	Conical
Self-Standing	Yes	Yes	Yes	Yes
Graduations	Yes	Yes	Yes	Yes
Graduation Range	5-50ml	5-50ml	5-50ml	5-50ml
Graduation Interval	2.5ml	2.5ml	2.5ml	2.5ml
Cap Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Cap Style	screw cap	screw cap	screw cap	screw cap
Cap Color	Yellow	Green	Yellow	Green
Max RCF	-	-	-	-
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Sterile	No	No	Yes	Yes
Dnase-free	No	No	Yes	Yes
Rnase-free	No	No	Yes	Yes
Outer Diameter	29mm	29mm	29mm	29mm
Length	117mm	117mm	117mm	117mm
Cap Diameter	33.3mm	33.3mm	33.3mm	33.3mm