



## **PRODUCT INFORMATION**

### **Quantitative Pipettes Fillers**

***Product information for PF330:***

**Description:**

- | Red Color
- | rubber
- | three valves
- | Natural rubber bulb allows safe, manual control of filling or discharging toxic liquids through a pipette by finger pressure (thus avoiding mouth contact). Guaranteed to be consistently flexible.
- | The maximum capacity is 50 ml.
- | PF330 can be attached to 5/10/25 ml pipette.



**Item type:**

Item No.	Type
PF330	Regular