



## CERTIFICATE OF ANALYSIS

---

Product	BES, sodium salt acid
Grade	High Purity
Product Code	BB0131
Formula	$C_6H_{14}NO_5SNa$
MW	235.2
CAS#	66992-27-6
Lot No	

---

---

Test Items	Specifications	Results
Appearance	White crystalline powder	
Purity	>99%	
DNase, RNase, Protease	None Detected	
Useful pH range	6.4~7.8	
pKa at 25°C	7.1 (free acid)	
Solubility	10g in 30ml water, clear colorless solution	
Loss on drying	<1%	
Heavy metals	<5ppm	

---

Storage: R.T.