



## CERTIFICATE OF ANALYSIS

---

Product	dNTP Mixture, 25mM
Grade	Biotech
Formula	N/A
M.W.	N/A
CAS#	N/A
Product Code	DD0057
Batch No	

---

---

Test Items	Specifications	Results
Apperance	Clear and colorless solution	
Dnase, Rnase	None detected	
pH	7.0 ±0.1	

---

Description: A solution that is 25mM dATP, 25mM dGTP, 25mM dCTP, 25mM dTTP.

Storage: Long term (infrequent use, 1-2 times per month): -70°C  
Daily/Weekly use: -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability.