



CERTIFICATE OF ANALYSIS

Product	Polypeptone
Product Code	G258
Lot	

Test Items	Specifications	Results
Appearance	Dry, fine, homogeneous, light to medium tan yellow power	
Total Nitrogen(TN)	≥ 10%	
Loss on Drying	≤ 6%	
Ash(Residue on ignition)	≤ 15%	
pH	6.5-7.5 at 25oC	
<u>Microbiological Analysis</u>		
Standard plate count	<5000 CFU/g	
Yeast and molds	<100 CFU/g	
Coliforms	Negative	
Salmonella	Negative	

Storage: R.T. Protect from moisture