



CERTIFICATE OF ANALYSIS

| | |
|--------------|---------------------------|
| Product | Okadaic acid, sodium salt |
| Grade | High Purity |
| Product Code | OL1199 |
| Formula | $C_{44}H_{67}NaO_{13}$ |
| MW | 826.98 |
| CAS# | 209266-80-8 |
| Lot No | |

| Test Items | Specifications | Results |
|------------|--|---------|
| Appearance | White, Crystalline solid or clear colorless solid film | |
| TLC purity | >98% | |
| Solubility | soluble in DMSO, methanol and ethanol | |

Storage: -20°C