



Bio Basic Inc.

A world leader in serving science

CERTIFICATE OF ANALYSIS

Product	PHOSPHATE BUFFERED SALINE(PBS)
Grade	BIOTECH GRADE
Product Code	PD0100
Lot No	

Test Items	Specifications	Results
Appearance	White granular powder	
Conductivity(1.0%, Water)	13500-17000umhos	
pH(1.0%, Water) @25C	7.3-7.5	
Solubility(1.0%, Water)	PASS	

1 x PBS solution contains 137 mM Sodium Chloride, 2.7mM Potassium Chloride and 10mM Phosphate Buffer Each pack makes 1 L of 10 concentrated PBS or 10 liters of 1 x solution

Storage: Store at room temperature. Protect from moisture & light.

QF 21 Rev 2015

20 Konrad Cres, Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com