



CERTIFICATE OF ANALYSIS

Product	LB BROTH TOP AGAR
Grade	N/A
Product Code	SD7004
Appearance:	Light beige, free-flowing powder
Lot No	

Test Items	Specification	Result
Conductivity (1:10) @25C	REPORT	
pH (6.4g/qs 200ml water)@25C	6.3-7.3	
Solubility (6.4g/qs 200ml water)	REPORT	
Supports Growth of Bacteria	PASS	

Component (g/L):

Tryptone	10.0g
Yeast extracts	5.0g
NaCl	10.0g
Agar	7.0g

Usage: 32.0g/L

Storage: Store at room temperature. Protect from moisture.