



CERTIFICATE OF ANALYSIS

Product	YPD-BROTH AGAR POWDER
Product Code	SD7023
Appearance:	Light beige, free-flowing powder
Lot No	

Test Items	Specifications	Results
Conductivity (1:10)@25C	REPORT	
pH (6.5g/ 100ml water) @25C	6.0-7.0	
Solubility (6.5g/ 100ml water)	PASS	
Supports Growth of Bacteria	PASS	

Component (g/L):

Yeast Extract	10.0g
Peptone	20.0g
D-glucose	20.0g
Agar	15.0g

Usage: 65.0g/L

Storage: Store at room temperature. Protect from moisture.