



CERTIFICATE OF ANALYSIS

Product	Uridine -5'-triphosphate UTP], 100 mM solution
Grade	High Purity
Product Code	UB1341
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	

Test Items	Specifications	Results
Appearance	Clear Colorless Solution	
Purity (HPLC)	≥99%	
Concentration	100mM ± 2%	
pH	6.8-7.2	

Description: UTP (Uridine -5'-triphosphate) is supplied as 100mM aqueous solution titrated to pH 7.0 with NaOH.

Application: For use in molecular biology applications in aRNA synthesis, siRNA synthesis, vitro transcription RNA amplification.

Storage: -20°C