



Product Information **EZ-LINE Strip Immunoplate, Megabinding**

Product information for SP49207:

Description

Plate contains twelve strips of eight wells each, for a total of 96 wells. Strips are placed in a white frame with alphanumeric labeling. Megabinding surface has a high affinity for proteins and other molecules containing both hydrophilic and hydrophobic regions, such as immunoglobulins.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- Labelled white frame
- High strength polystyrene construction
- Megabinding surface



Item No.	Description	Package
SP49207	EZ-LINE Strip Immunoplate, 8 well/strip, Megabinding, 85.4x127.6mm sterile, 10 per pack, 100 per case	100/case

Technical Specifications

Item No.	SP49207
Quantity per case	100/case
Packaging	10 Pack
Volume	350uL
Growth Area	0.33cm ²
Color	Clear Strips, White Frame
Plastic materials	Polystyrene
Sterile	Yes
Surface Type	Megabinding

Storage

Room temperature.