



CERTIFICATE OF ANALYSIS

Product	LB Agar Powder, Miller
Product Code	SD7003(S519)
CAS#	N/A
Lot#	

Test Items	Specifications	Results
Appearance	Off-White to Light Brown Powder	
Solubility (3.7g/ 100ml, H ₂ O)	Clear to slightly hazy, Faint yellow to yellow solution	

Storage: 18 ~25°C. Protect from moisture

Composition (g/L)

Trypton	10g
Yeast Extract	5g
Sodium Chloride	10g
Agar	15g

Usage: 40g/L