



Product Information

MCE Membrane

Product information for HP132-HP135:

Description

Membranes (code# HP132 to HP135) are for aqueous solution. They are made of mixed fiber (MCE). They are suitable for filtration in various aqueous solutions. The membrane is not resistant to organic liquid and strong acid and strong alkali solution. The membrane can withstand 121°C for 20 minutes.



| Item No. | Membrane | Diameter | Pore Size | Package |
|-----------------|-----------------|-----------------|------------------|----------------|
| HP132 | MCE | 25mm | 0.45um | 100/box |
| HP133 | MCE | 50mm | 0.45um | 100/box |
| HP134 | MCE | 50mm | 0.22um | 100/box |
| HP135 | MCE | 25mm | 0.80um | 100/box |

Storage

Room temperature.