



PRODUCT INFORMATION

PSB Tablet (PD0435)

Product name: PBS (Phosphate buffered saline),

Catalog Number: PD0435

Application:

Phosphate Buffered Saline tablets have been specifically developed for use in immunoassay procedures.

Reconstitution:

One tablet dissolved in 100 ml of deionized water yields 1 X PBS buffer. The final 1X solution contains: 2.7 mM potassium chloride, 137 mM sodium chloride and 10 mM Phosphate Buffer, pH: 7.3-7.5 at 25 °C (1 Tablet in 100ml water).

Storage Temperature: Room temperature.