



CERTIFICATE OF ANALYSIS

Product	PBS (Phosphate buffered saline), 100 ml tablets
Grade	Biotech
Product Code	PD0435(D0435)
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	

Test Items	Specifications	Results
Appearance	White Round Tablet	
pH (1 tablet/100mL Water at 25°C)	7.3~7.5	
Usage	Each tablet prepares 100 ml of a 1X solution.	

1X PBS Solution contains 137 mM Sodium Chloride, 2.7 mM Potassium Chloride and 10 mM Phosphate Buffer.

Storage: R.T.